

**bs-15244R**

**[ Primary Antibody ]**

## C6orf226 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |                                                        |                                                                                                                                                                                                                                                                                                                                                              |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>GeneID:</b> 441150</p> <p><b>Target:</b> C6orf226</p> <p><b>Immunogen:</b> KLH conjugated synthetic peptide derived from human C6orf226 : 1-80/101.</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.<br/>Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> | <p><b>Isotype:</b> IgG</p> <p><b>SWISS:</b> Q5I0X4</p> | <p><b>Applications:</b> <b>WB</b> (1:500-2000)<br/><b>IHC-P</b> (1:100-500)<br/><b>IHC-F</b> (1:100-500)<br/><b>IF</b> (1:100-500)<br/><b>ICC/IF</b> (1:100-500)<br/><b>ELISA</b> (1:5000-10000)</p> <p><b>Reactivity:</b> (predicted: Human, Rat)</p> <p><b>Predicted MW.:</b> 11 kDa</p> <p><b>Subcellular Location:</b> Extracellular matrix ,Nucleus</p> |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|