
Kazrin Rabbit pAb

Catalog Number: bs-9707R

Target Protein: Kazrin

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:50-200), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Pig, Cow, Dog)

Predicted MW: 86 kDa

Entrez Gene: 23254

Swiss Prot: Q674X7

Source: KLH conjugated synthetic peptide derived from human Kazrin: 451-550/775.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Component of the cornified envelope of keratinocytes. May be involved in the interplay between adherens junctions and desmosomes. The function in the nucleus is not known.