## bs-9707R

## [ Primary Antibody ]

www.bioss.com.cn

sales@bioss.com.cn techsupport@bioss.com.cn

400-901-9800

DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Kazrin Rabbit pAb

GenelD: 23254 **SWISS:** Q674X7

Target: Kazrin

**Immunogen:** KLH conjugated synthetic peptide derived from human Kazrin:

451-550/775.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Background: Component of the cornified envelope of keratinocytes. May be

involved in the interplay between adherens junctions and desmosomes. The function in the nucleus is not known.

Applications: WB (1:500-2000)

**IHC-P** (1:100-500) **IHC-F** (1:100-500) **IF** (1:100-500)

**ELISA** (1:5000-10000)

Reactivity: (predicted: Human, Mouse,

Rat, Pig, Cow, Dog)

Predicted 86 kDa

**Subcellular** Cytoplasm