

## Cardiac Troponin I-C (cTnI-C) Mouse mAb

|                      |                                                                                 |
|----------------------|---------------------------------------------------------------------------------|
| Cat. No:             | V3501                                                                           |
| Name:                | Cardiac Troponin I-C (cTnI-C) Mouse mAb                                         |
| Storage:             | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type:          | Monoclonal                                                                      |
| Clone No:            | 5D12K                                                                           |
| Host/Source:         | Mouse                                                                           |
| Subtype:             | IgG1                                                                            |
| Reactivity:          | Human                                                                           |
| Format:              | Liquid                                                                          |
| Concentration:       | >1mg/ml                                                                         |
| Buffer:              | 0.01M PBS (pH7.4) with 0.02% Proclin300.                                        |
| Applications:        | ELISA                                                                           |
| Antigen:             | Recombinant human Cardiac Troponin I-C ( <a href="#">bs-41212P</a> )            |
| Pairing Information: |                                                                                 |

| ELISA   |               | Detection |       |       |               |               |
|---------|---------------|-----------|-------|-------|---------------|---------------|
|         |               | V3201     | V3401 | V3501 | V3201 + V3401 | V3401 + V3501 |
| Capture | V3201         | -         | -     | +++   | -             | -             |
|         | V3401         | +++       | -     | -     | -             | -             |
|         | V3501         | -         | -     | -     | -             | -             |
|         | V3201 + V3401 | -         | -     | +++   | -             | -             |
|         | V3401 + V3501 | +++       | -     | -     | -             | -             |