
MTMR10 Rabbit pAb

Catalog Number: bs-17882R

Target Protein: MTMR10

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Pig, Sheep, Dog, Horse)

Predicted MW: 88 kDa

Entrez Gene: 54893

Swiss Prot: Q9NXD2

Source: KLH conjugated synthetic peptide derived from human MTMR10: 261-360/777.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Probable pseudophosphatase. Contains a Glu residue instead of a conserved Cys residue in the dsPTPase catalytic loop which renders it catalytically inactive as a phosphatase.