

bs-0991R

[Primary Antibody]

HPV18 E6 protein Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: HPV18 E6 protein

Immunogen: KLH conjugated synthetic peptide derived from HPV18 E6: 75-158/158.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ELISA (1:5000-10000)

Reactivity: (predicted: Human papillomavirus type 18)

Predicted MW.: 11 kDa