## bs-0991R

## [ Primary Antibody ]

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

## **HPV18 E6 protein Rabbit pAb**

- DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Target: HPV18 E6 protein

**Immunogen:** KLH conjugated synthetic peptide derived from HPV18 E6: 75-158/158.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

**Applications: IHC-P** (1:100-500)

**IHC-F** (1:100-500) **IF** (1:100-500) **ELISA** (1:5000-10000)

Reactivity: (predicted: HPV18)

Predicted MW.: 11 kDa