

**bs-7002R**

**[ Primary Antibody ]**

## PAGE4 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>GeneID:</b> 9506</p> <p><b>Target:</b> PAGE4</p> <p><b>Immunogen:</b> KLH conjugated synthetic peptide derived from human PAGE4: 45-102/102.</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p><b>Background:</b> PAGE4 belongs to the GAGE family of genes which are expressed in a variety of tumors and in some fetal and reproductive tissues. PAGE4 is strongly expressed in prostate and prostate cancer, but is also expressed in other male and female reproductive tissues including testis, fallopian tube, uterus, and placenta, as well as in testicular cancer and uterine cancer.</p>	<p><b>Isotype:</b> IgG</p> <p><b>Applications:</b> <b>ELISA</b> (1:5000-10000)</p> <p><b>Reactivity:</b> (predicted: Human)</p> <p><b>Predicted MW.:</b> 11 kDa</p> <p><b>Subcellular Location:</b> Cytoplasm ,Nucleus</p>
---	--