

bs-0925R

[Primary Antibody]

Vibrio cholerae O1 Serotype Inaba Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Applications: ELISA (1:5000-10000)

Clonality: Polyclonal

Target: Vibrio cholerae O1 Serotype Inaba

Reactivity: (predicted: Vibrio
cholerae O1 Serotype
Inaba)

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.