

**bs-21263R**

**[ Primary Antibody ]**

## SLC35E2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<b>Host:</b> Rabbit <b>Clonality:</b> Polyclonal <b>GeneID:</b> 55508 <b>Target:</b> SLC35E2 <b>Purification:</b> affinity purified by Protein A <b>Concentration:</b> 1mg/ml <b>Storage:</b> Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.	<b>Isotype:</b> IgG <b>SWISS:</b> P0CK97	<b>Applications:</b> WB (1:500-2000) <b>Reactivity:</b> Mouse (predicted: Human, Rat, Pig, Sheep, Cow, Dog, Horse) <b>Predicted MW.:</b> 29 kDa <b>Subcellular Location:</b> Cell membrane
----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------