

---

## SLC35E2 Rabbit pAb

Catalog Number: bs-21263R

Target Protein: SLC35E2

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000)

Reactivity: Mouse (predicted:Human, Rat, Pig, Sheep, Cow, Dog, Horse)

Predicted MW: 29 kDa

Entrez Gene: 55508

Swiss Prot: P0CK97

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.