
ZDHHC1/ZNF377 Rabbit pAb

Catalog Number: bs-18471R

Target Protein: ZDHHC1/ZNF377

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500)

Reactivity: (predicted:Human, Mouse, Rat)

Predicted MW: 55 kDa

Subcellular: Cell membrane

Locations:

Entrez Gene: 29800

Swiss Prot: Q8WTX9

Source: KLH conjugated synthetic peptide derived from human ZDHHC1/ZNF377: 1-100/485.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.