

**bs-18471R**

**[ Primary Antibody ]**

## **ZDHHC1/ZNF377 Rabbit pAb**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### **— DATASHEET —**

<b>Host:</b> Rabbit	<b>Isotype:</b> IgG	<b>Applications:</b> <b>WB</b> (1:500-2000)
<b>Clonality:</b> Polyclonal		<b>IHC-P</b> (1:100-500)
<b>GeneID:</b> 29800	<b>SWISS:</b> Q8WTX9	<b>IHC-F</b> (1:100-500)
<b>Target:</b> ZDHHC1/ZNF377		<b>IF</b> (1:100-500)
<b>Immunogen:</b> KLH conjugated synthetic peptide derived from human ZDHHC1/ZNF377: 1-100/485.		<b>ICC/IF</b> (1:100-500)
<b>Purification:</b> affinity purified by Protein A		<b>Reactivity:</b> (predicted: Human, Mouse, Rat)
<b>Concentration:</b> 1mg/ml		
<b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		<b>Predicted MW.:</b> 55 kDa
		<b>Subcellular Location:</b> Cell membrane