

bs-16961R

[Primary Antibody]

NS5ATP9 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 9768</p> <p>Target: NS5ATP9</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human NS5ATP9: 21-80/111.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: NS5ATP9/PAF/p15(PAF) is involved in protection of cells from UV-induced cell death and is also a Hepatitis C NS5A up-regulation gene which may play a role in the pathogenesis of HCV-associated hepatocellular carcinoma.</p>	<p>Isotype: IgG</p> <p>SWISS: Q15004</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Pig, Sheep, Cow, Dog, Horse)</p> <p>Predicted MW.: 12 kDa</p> <p>Subcellular Location: Nucleus</p>
---	--	---