

Human S100B Mouse mAb

Cat. No: V2501
Name: Human S100B Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 18F6
Host/Source: Mouse
Subtype: IgG1
Reactivity: (predicted: Human)
Format: Liquid
Concentration: >1mg/ml
Buffer: 0.01M PBS (pH7.4) with 0.02% Proclin300.
Applications: ELISA, ICA
Antigen: Recombinant human S100B protein ([bs-41191P](#))
Pairing Information:

ELISA		Detection			
		V2501	V2503	V2504	V2505
Capture	V2501	-	-	-	-
	V2503	+++	-	-	+
	V2504	-	+	-	-
	V2505	+	+	-	-
ICA		Detection			
		V2501	V2503	V2504	V2505
Capture	V2501	-	-	-	-
	V2503	+	-	+++	-
	V2504	-	-	-	-
	V2505	-	-	-	-