

Human SAA1 Mouse mAb

Cat. No: V0201
Name: Human SAA1 Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 3A8
Host/Source: Mouse
Subtype: IgG2b
Reactivity: (predicted: Human)
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA, ICA
Antigen: Recombinant Human Serum Amyloid A protein ([bs-41127P](#))
Pairing Information:

ELISA		Detection	
		V0201	V0203
Capture	V0201	-	+++
	V0203	-	-
ICA		Detection	
		V0201	V0203
Capture	V0201	-	-
	V0203	+++	-