
SAA1 Mouse mAb

Cat. No:	V0202
Name:	SAA1 Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted: Human)
Format:	Liquid
Concentration:	1mg/ml
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA
Antigen:	Recombinant Human Serum Amyloid A protein (bs-41127P)