

Human SAA1 Mouse mAb

Cat. No:	V0203
Name:	Human SAA1 Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	8A3
Host/Source:	Mouse
Subtype:	IgM
Reactivity:	(predicted: Human)
Format:	Liquid
Concentration:	1mg/ml
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA, ICA
Antigen:	Recombinant Human Serum Amyloid A protein (bs-41127P)
Pairing Information:	

ELISA		Detection	
		V0201	V0203
Capture	V0201	-	+++
	V0203	-	-
ICA		Detection	
		V0201	V0203
Capture	V0201	-	-
	V0203	+++	-