

**bsm-0910M**

**[ Primary Antibody ]**

## Batroxobin(5H4) Mouse mAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

**Host:** Mouse

**Isotype:** IgG

**Clonality:** Monoclonal

**CloneNo.:** 5H4

**Target:** Batroxobin(5H4)

**Purification:** affinity purified by Protein G

**Concentration:** 1mg/ml

**Storage:** Size : 50ul/100ul/200ul  
0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.  
Size : 200ug (PBS only)  
0.01M PBS  
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Background:** No data available.

**Applications:** IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ELISA (1:5000-10000)

**Reactivity:** (predicted: Agkistrodon halys)

**Predicted MW.:** 25 kDa