

bs-16969R

[Primary Antibody]

KIAA0323 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) ELISA (1:5000-10000)
Clonality: Polyclonal		
GeneID: 23351	SWISS: O15037	Reactivity: (predicted: Human, Mouse, Rat)
Target: KIAA0323		
Immunogen: KLH conjugated synthetic peptide derived from human KIAA0323: 461-550/678.		Predicted MW.: 75 kDa
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: KIAA0323 belongs to the N4BP1 family. Its function is not known.		