

bs-6282R

[Primary Antibody]

SERPINB4 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Dog, Horse) Predicted MW.: 45 kDa Subcellular Location: Cytoplasm
Clonality: Polyclonal		
GeneID: 6318	SWISS: P48594	
Target: SERPINB4		
Immunogen: KLH conjugated synthetic peptide derived from human SERPINB4: 121-220/390.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: May act as a protease inhibitor to modulate the host immune response against tumor cells.		

— SELECTED CITATIONS —

- **[IF=6.057]** Limin Zhou. et al. Copper doped terbium metal organic framework as emitter for sensitive electrochemiluminescence detection of CYFRA 21-1. Talanta. 2022 Feb;238:123047 Other ;. 34801904