
SERPINB4 Rabbit pAb

Catalog Number: bs-6282R
Target Protein: SERPINB4
Concentration: 1mg/ml
Form: Liquid
Host: Rabbit
Clonality: Polyclonal
Isotype: IgG
Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ELISA (1:5000-10000)
Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Dog, Horse)
Predicted MW: 45 kDa
Entrez Gene: 6318
Swiss Prot: P48594
Source: KLH conjugated synthetic peptide derived from human SERPINB4: 121-220/390.
Purification: affinity purified by Protein A
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Background: May act as a protease inhibitor to modulate the host immune response against tumor cells.

PRODUCT SPECIFIC PUBLICATIONS

[IF=6.057] Limin Zhou. et al. Copper doped terbium metal organic framework as emitter for sensitive electrochemiluminescence detection of CYFRA 21-1. Talanta. 2022 Feb;238:123047 Other ; . 34801904