
IBRDC1 Rabbit pAb

Catalog Number: bs-15531R

Target Protein: IBRDC1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500),
ELISA (1:5000-10000)

Reactivity: (predicted:Human, Rat, Sheep, Cow, Dog, Horse)

Predicted MW: 32 kDa

Subcellular Cell membrane

Locations:

Entrez Gene: 154214

Swiss Prot: Q8TC41

Source: KLH conjugated synthetic peptide derived from human IBRDC1: /51-150/275.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

PRODUCT SPECIFIC PUBLICATIONS

[IF=3.2] Wanqing Wang. et al. STL Inhibited Angiogenesis of DPSCs Through Depressing Mitochondrial Respiration by Enhancing RNF217.

Advanced Biology. 2024 Jun;;2400042 WB ; Human . 38880848