bs-17105R

[Primary Antibody]

LARP6 Rabbit pAb

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

GenelD: 55323 SWISS: Q9BRS8

Target: LARP6

Immunogen: KLH conjugated synthetic peptide derived from human LARP6:

151-250/419.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500) **IF** (1:100-500) ICC/IF (1:100-500)

Reactivity: (predicted: Human, Mouse,

Rat, Pig, Sheep, Cow,

Horse)

Predicted 55 kDa

Subcellular Location: Cytoplasm ,Nucleus

- SELECTED CITATIONS -

• [IF=3.9] Kong Lingfan. et al. IGF-1 regulates LARP6-mediated collagen metabolism in vaginal fibroblasts of POP patients via the PI3K/AKT pathway. SCI REP-UK. 2025 Jul;15(1):1-20 IHC; Human. 40593278