

bs-6754R

[Primary Antibody]

YY1AP1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 55249</p> <p>Target: YY1AP1</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human YY1AP1/HCCA2: 301-400/796.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: Enhances transcription activation by YY1. Tissue specificity: Ubiquitous. Detected in small intestine, skeletal muscle, lung, pancreas, brain, stomach, spleen, colon and heart. Detected at very low levels in healthy liver. Highly expressed in most liver carcinomas.</p>	<p>Isotype: IgG</p> <p>SWISS: Q9H869</p>	<p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human)</p> <p>Predicted MW.: 88 kDa</p> <p>Subcellular Location: Cytoplasm ,Nucleus</p>
---	--	--