

---

## YY1AP1 Rabbit pAb

Catalog Number: bs-6754R

Target Protein: YY1AP1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: **WB** (1:500-2000), **IHC-P** (1:100-500), **IHC-F** (1:100-500), **IF** (1:100-500), **ELISA** (1:5000-10000)

Reactivity: (predicted:&nbsp;Human)

Predicted MW: 88 kDa

Subcellular Cytoplasm ,Nucleus

Locations:

Entrez Gene: 55249

Swiss Prot: Q9H869

Source: KLH conjugated synthetic peptide derived from human YY1AP1/HCCA2: 301-400/796.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Background:** Enhances transcription activation by YY1. Tissue specificity:Ubiquitous. Detected in small intestine, skeletal muscle, lung, pancreas, brain, stomach, spleen, colon and heart. Detected at very low levels in healthy liver. Highly expressed in most liver carcinomas.