

---

## NRIP1 Rabbit pAb

Catalog Number: bs-8739R

Target Protein: NRIP1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Pig, Sheep, Cow, Horse)

Predicted MW: 140 kDa

Entrez Gene: 8204

Swiss Prot: P48552

Source: KLH conjugated synthetic peptide derived from human NRIP1: 1001-1100/1158.

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.