

bs-5836R

[Primary Antibody]

EMP3 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: ELISA (1:5000-10000)
Clonality: Polyclonal		Reactivity: Human (predicted: Mouse, Rat, Pig, Sheep, Cow, Horse)
GeneID: 2014	SWISS: P54852	Predicted MW.: 18 kDa
Target: EMP3		Subcellular Location: Cell membrane
Immunogen: KLH conjugated synthetic peptide derived from human EMP3: 101-163/163.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: EMP3 is an integral membrane protein that may be involved in cell proliferation and cell-cell interactions.		