

bs-11772R

[Primary Antibody]

Tetanus Toxin heavy chain Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: Tetanus Toxin heavy chain

Immunogen: KLH conjugated synthetic peptide derived from Clostridium tetani E88 Tetanus toxin heavy chain: 801-900/1315.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Applications: IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ICC/IF (1:100-500)

ELISA (1:5000-10000)

Reactivity: (predicted: Clostridium tetani E88)

Predicted MW.: 94 kDa