

bs-5826R

[Primary Antibody]

PCDHAC1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

| | | |
|--|--|--|
| <p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 56135</p> <p>Target: PCDHAC1</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human PCDHAC1: 551-650/963. < Extracellular ></p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: Potential calcium-dependent cell-adhesion protein. May be involved in the establishment and maintenance of specific neuronal connections in the brain.</p> | <p>Isotype: IgG</p> <p>SWISS: Q9H158</p> | <p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat)</p> <p>Predicted MW.: 102 kDa</p> <p>Subcellular Location: Cell membrane</p> |
|--|--|--|