## bs-17182R

## [ Primary Antibody ]

## Bioss

## TRAF3IP1 Rabbit pAb

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

GeneID: 26146 SWISS: Q8TDR0

Target: TRAF3IP1

**Immunogen:** KLH conjugated synthetic peptide derived from human TRAF3IP1:

501-600/691.

Purification: affinity purified by Protein A

Concentration: 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

**Applications: IHC-P** (1:100-500)

IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)

Reactivity: (predicted: Human, Mouse,

Rat, Rabbit, Pig, Sheep, Cow, Zebrafish, Chicken, Dog, GuineaPig, Horse)

Predicted MW.: 79 kDa

Subcellular Location: Cytoplasm