

bs-4232R

[Primary Antibody]

ATG16L2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

| | | |
|--|----------------------|---|
| Host: Rabbit | Isotype: IgG | Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse) Predicted MW.: 69 kDa Subcellular Location: Cytoplasm |
| Clonality: Polyclonal | | |
| GeneID: 89849 | SWISS: Q8NAA4 | |
| Target: ATG16L2 | | |
| Immunogen: KLH conjugated synthetic peptide derived from human ATG16L2: 191-290/619. | | |
| Purification: affinity purified by Protein A | | |
| Concentration: 1mg/ml | | |
| Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. | | |
| Background: May play a role in autophagy. | | |