

bs-5831R

[Primary Antibody]

ANKHD1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Pig, Sheep, Cow, Dog, Horse) Predicted MW.: 269 kDa
Clonality: Polyclonal		
GeneID: 54882	SWISS: Q8IWZ3	
Target: ANKHD1		
Immunogen: KLH conjugated synthetic peptide derived from human ANKHD1: 1901-2000/2542.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: May play a role as a scaffolding protein that may be associated with the abnormal phenotype of leukemia cells. Isoform 2 may possess an antiapoptotic effect and protect cells during normal cell survival through its regulation of caspases.		

— SELECTED CITATIONS —

- **[IF=6.083]** Sun, Wenhui. et al. Alternative ANKHD1 transcript promotes proliferation and inhibits migration in uterine corpus endometrial carcinoma. NPJ GENOM MED. 2022 Sep;7(1):1-12 IF ;Human. 36171217