

bs-6777R

[Primary Antibody]

HIPK3 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

| | | |
|--|--|---|
| <p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 10114</p> <p>Target: HIPK3</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human HIPK3: 401-500/1215.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: Seems to negatively regulate apoptosis by promoting FADD phosphorylation. Enhances androgen receptor-mediated transcription. May act as a transcriptional corepressor for NK homeodomain transcription factors.</p> | <p>Isotype: IgG</p> <p>SWISS: Q9H422</p> | <p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Sheep, Cow, Chicken, Dog, Horse)</p> <p>Predicted MW.: 134 kDa</p> <p>Subcellular Location: Cytoplasm ,Nucleus</p> |
|--|--|---|