

bs-0312R-Gold

[Secondary Antibodies]



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

Rabbit Anti-Avidin, Gold conjugated

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: Rabbit Anti-Avidin

Purification: affinity purified by Protein A

Concentration: 0.4 mg/ml

Storage: 20mM TBS (pH8.0) with 1% BSA and 0.02% Proclin300.
Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.

Background: The biological function of avidin is not known. It is a homotetrameric glycoprotein found in chicken egg white, and forms a strong non covalent specific complex with biotin (one molecule of biotin per subunit of avidin).

Applications: **ICA** (1:20-200)
IEM (1:20-200)

Reactivity: Avidin