
Rabbit Anti-Avidin, Gold conjugated

Catalog Number: bs-0312R-Gold

Target Protein: Rabbit Anti-Avidin

Concentration: 0.4 mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Reactivity: Avidin

Purification: affinity purified by Protein A

Storage: 20mM TBS (pH8.0) with 1% BSA and 0.02% Proclin300.

Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.

Background: The biological function of avidin is not known. It is a homotetrameric glycoprotein found in chicken egg white, and forms a strong non covalent specific complex with biotin (one molecule of biotin per subunit of avidin).