
HBLD1 Rabbit pAb

Catalog Number: bs-15416R

Target Protein: HBLD1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Cow, Dog, Horse)

Predicted MW: 16 kDa

Subcellular Cytoplasm

Locations:

Entrez Gene: 122961

Swiss Prot: Q86U28

Source: KLH conjugated synthetic peptide derived from human HBLD1: 61-154/154.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.