

bs-15416R

[Primary Antibody]

HBLD1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: ELISA (1:5000-10000)
Clonality: Polyclonal		Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Cow, Dog, Horse)
GeneID: 122961	SWISS: Q86U28	Predicted MW.: 16 kDa
Target: HBLD1		Subcellular Location: Cytoplasm
Immunogen: KLH conjugated synthetic peptide derived from human HBLD1: 61-154/154.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		