

bs-0903R

[Primary Antibody]

Naltrexone Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: Naltrexone

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Naltrexone is a drug prescribed to assist people in maintaining abstinence after they have detoxified or withdrawn from opioids (eg heroin) or alcohol. Naltrexone is classed as an opioid agonist. It works by blocking the opioid receptors in the brain so that opiates have no effect on the brain.

Applications: ELISA (1:5000-10000)

Reactivity: (predicted: Naltrexone)