

bs-0331R

[Primary Antibody]

Hepatitis B Virus Surface Antigen Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Polyclonal

Target: Hepatitis B Virus Surface Antigen

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background:

Applications: IHC-P (1:100-500)

IHC-F (1:100-500)

IF (1:100-500)

ICC/IF (1:100-500)

Reactivity: (predicted: Hepatitis B virus (HBVE1))

Predicted MW.: 44 kDa

Subcellular Location: Cell membrane ,Cytoplasm