

bs-14518R

[Primary Antibody]

EFHA2 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 286097</p> <p>Target: EFHA2</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: MICU3 (Mitochondrial Calcium Uptake Family Member 3) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include calcium ion binding. An important paralog of this gene is MICU2.</p>	<p>Isotype: IgG</p> <p>SWISS: Q86XE3</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Sheep, Cow, Dog, Horse)</p> <p>Predicted MW.: 61 kDa</p> <p>Subcellular Location: Cytoplasm</p>
---	--	--