
Procalcitonin (PCT) Mouse mAb

Cat. No:	V0305
Name:	Procalcitonin (PCT) Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3E3L
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted: Human)
Format:	Liquid
Concentration:	≥1mg/ml
Buffer:	0.01M PBS (pH7.4) with 0.02% Proclin300.
Applications:	ELISA
Antigen:	Recombinant human Procalcitonin (bs-41144P)