bs-0343R

[Primary Antibody]

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn

- DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Target: large S protein

Immunogen: KLH conjugated synthetic peptide derived from HBSAG_HBVC5:

201-281/281.

large S protein Rabbit pAb

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Background:

Applications: WB (1:500-2000)

400-901-9800

ELISA (1:5000-10000)

Reactivity: (predicted: Hepatitis B

Predicted MW.: 31/44 kDa

Subcellular Location: Cell membrane, Cytoplasm

- SELECTED CITATIONS -

• [IF=1.93] Yuan, Xiaocong. "A label-free approach to kinetic analysis and high multiplex detection of targeted drugs with phase surface plasmon resonance imaging. "Analytical Methods (2015). Other; 10.1039/C4AY02174F