

bs-0343R

[Primary Antibody]

large S protein Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) ELISA (1:5000-10000)
Clonality: Polyclonal		Reactivity: (predicted: Hepatitis B virus)
Target: large S protein		Predicted MW.: 31/44 kDa
Immunogen: KLH conjugated synthetic peptide derived from HBSAG_HBVC5: 201-281/281.		Subcellular Location: Cell membrane ,Cytoplasm
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background:		

— SELECTED CITATIONS —

- **[IF=1.93]** Yuan, Xiacong. "A label-free approach to kinetic analysis and high multiplex detection of targeted drugs with phase surface plasmon resonance imaging."Analytical Methods (2015). Other ;="". 10.1039/C4AY02174F