
OSGIN2 Rabbit pAb

Catalog Number: bs-23505R

Target Protein: OSGIN2

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500)

Reactivity: (predicted:Human, Mouse, Rat)

Predicted MW: 57 kDa

Entrez Gene: 734

Swiss Prot: Q9Y236

Source: KLH conjugated synthetic peptide derived from human OSGIN2 : 131-230/505.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: OSGIN2 (Oxidative stress induced growth inhibitor family member 2) may be involved in meiosis or the maturation of germ cells.