## bs-12019R

## [ Primary Antibody ]

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

## **GPR1** Rabbit pAb

- DATASHEET -

Isotype: IgG Host: Rabbit

Clonality: Polyclonal

Target: GPR1

**Immunogen:** KLH conjugated synthetic peptide derived from human GPR1: 151-250/355. < Extracellular >

Purification: affinity purified by Protein A

Concentration: 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

**Applications: ELISA** (1:5000-10000)

Reactivity: (predicted: Human,

Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog,

Horse)

Predicted MW.: 41 kDa