

bsm-41483M**[Primary Antibody]****Bioss**
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800


SARS-CoV-2 (2019-nCoV) Nucleocapsid Mouse mAb**DATASHEET****Host:** Mouse**Isotype:** IgG**Clonality:** Monoclonal**CloneNo.:** 9D2-B2**Target:** SARS-CoV-2 (2019-nCoV) Nucleocapsid**Immunogen:** Recombinant SARS-CoV-2 Nucleocapsid Protein: 1-419/419.**Purification:** affinity purified by Protein A**Concentration:** 1mg/ml**Storage:** Size : 100ug

0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Size : 200ug (PBS only)

0.01M PBS

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: The SARS-CoV-2 spike (S) protein is the target of vaccine design efforts to end the COVID-19 pandemic. Despite a low mutation rate, isolates with the D614G substitution in the S protein appeared early during the pandemic, and are now the dominant form worldwide. Here, we analyze the D614G mutation in the context of a soluble S ectodomain construct.**Applications:** WB (1:500-2000)**Reactivity:** SARS-CoV-2**Predicted MW.:** 46 kDa**VALIDATION IMAGES**

Sample: Lane 1: Recombinant SARS-CoV-2 N protein (WT) (His Tag) (bs-41408P) Lane 2: Recombinant SARS-CoV-2 N protein (Q9H, P67S, P80R, P151L, S183Y) (His Tag) (bs-41451P) Lane 3: Recombinant SARS-CoV-2 N protein (D3L, P13T, D103Y, D128Y, H145Y, R203K, G204R, T205I, S235F) (His Tag) (bs-41452P) Lane 4: Recombinant SARS-CoV-2 N protein (Del204, Del215) (His Tag) (bs-41491P) Lane 5: Recombinant SARS-CoV-2 N protein (R203M, D377Y) (His Tag) (bs-41492P) Lane 6: Recombinant SARS-CoV-2 (Omicron, B.1.1.529) N protein (P13L, E31del, R32del, S33del, R203K, G204R) (N-His Tag) (bs-41494P) Primary: Anti-SARS-CoV-2(2019-nCoV)Nucleocapsid (bsm-41483M) at 1/1000 dilution Secondary: IRDye800CW Goat Anti-Mouse IgG at 1/20000 dilution Predicted band size: 46 kDa Observed band size: 50 kDa