bs-15384R

[Primary Antibody]

www.bioss.com.cn sales@bioss.com.cn

techsupport@bioss.com.cn 400-901-9800

DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

GPR158 Rabbit pAb

GenelD: 57512 **SWISS:** Q5T848

Target: GPR158

Immunogen: KLH conjugated synthetic peptide derived from human GPR158:

201-300/1215. < Extracellular >

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Applications: WB (1:500-2000)

IHC-P (1:100-500) **IHC-F** (1:100-500) **IF** (1:100-500) ICC/IF (1:100-500) **ELISA** (1:5000-10000)

Reactivity: (predicted: Human, Mouse,

Rat, Pig, Sheep, Cow,

Horse)

Predicted MW.: 133 kDa

Subcellular Location: Cell membrane

— SELECTED CITATIONS -

• [IF=6.2] Mu Lin. et al. Trilobatin, a Naturally Occurring GPR158 Ligand, Alleviates Depressive-like Behavior by Promoting Mitophagy. J AGR FOOD CHEM. 2025;73(9):5163-5179 IF; Mouse. 39962827