

## Digoxin Mouse mAb

Cat. No: V5101  
Name: Digoxin Mouse mAb  
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  
Clone type: Monoclonal  
Clone No: 6H7  
Host/Source: Mouse  
Subtype: IgG1  
Format: Liquid  
Concentration: >1mg/ml  
Buffer: 0.01M PBS (pH7.4) with 0.02% Proclin300.  
Applications: ELISA, ICA  
Antigen: KLH conjugated Digoxin  
Pairing Information:

ELISA		Detection	
		V5101	bs-0356PB
Capture	V5101	-	-
	bs-0356PB	+++	-
ICA		Detection	
		V5101	bs-0356PB
Capture	V5101	-	-
	bs-0356PB	+++	-