
D-Dimer Mouse mAb

Cat. No:	V1105
Name:	D-Dimer Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	2F7
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	Human
Format:	Liquid
Concentration:	≥0.6 mg/mg
Buffer:	potassium phosphate buffer and sodium chloride. No preservatives are added.
Applications:	ICA
Antigen:	Human fibrinogen cleavage product D-Dimer (bs-41172P)
Pairing Information:	
V1104-V1103&V1105	